

Clappac klapblisters

| Product | | | Afmeting | | | | | Foliemateriaal | | Transportverpakking | | | | |
|------------|-------------------|---------------|------------|-----|---------|-----|----|----------------|------|---------------------|--------------|------------|---------|-----------------|
| Artikel | Omschrijving | Type | Buitenmaat | | Vulmaat | | | Dikte MICRON | Type | Per doos | Aantal dozen | Per pallet | Gewicht | Gewichtseenheid |
| | | | B | H | B | H | D | | | | | | | |
| B042262-20 | Clamshell blister | Klapblister 1 | 192 | 82 | 157 | 62 | 20 | 370 | PET | 500 | 12 | 6.000 | 18,00 | Gram (g) |
| B042262-40 | Clamshell blister | Klapblister 1 | 192 | 82 | 157 | 62 | 40 | 370 | PET | 500 | 10 | 4.800 | 18,00 | Gram (g) |
| B112236-10 | Clamshell blister | Klapblister 1 | 92 | 57 | 66 | 48 | 10 | 310 | PET | 1400 | 12 | 16.800 | 10,00 | Gram (g) |
| B112236-15 | Clamshell blister | Klapblister 1 | 95 | 57 | 66 | 48 | 15 | 310 | PET | 1400 | 12 | 16.800 | 10,00 | Gram (g) |
| B112236-20 | Clamshell blister | Klapblister 1 | 92 | 57 | 66 | 48 | 20 | 330 | PET | 1400 | 12 | 16.800 | 10,00 | Gram (g) |
| B112236-30 | Clamshell blister | Klapblister 1 | 92 | 57 | 66 | 48 | 30 | 350 | PET | 1400 | 12 | 16.800 | 10,00 | Gram (g) |
| B112236-40 | Clamshell blister | Klapblister 1 | 92 | 57 | 66 | 48 | 40 | 420 | PET | 1400 | 12 | 16.800 | 10,00 | Gram (g) |
| B112408-30 | Clamshell blister | Klapblister 1 | 118 | 132 | 84 | 120 | 30 | 480 | PET | 400 | 12 | 4.800 | 14,00 | Gram (g) |
| B112408-40 | Clamshell blister | Klapblister 1 | 118 | 132 | 84 | 120 | 40 | 480 | PET | 400 | 12 | 4.800 | 14,00 | Gram (g) |
| B112666-25 | Clamshell blister | Klapblister 1 | 90 | 55 | 50 | 40 | 20 | 330 | PET | 1400 | 12 | 16.800 | 13,00 | Gram (g) |
| B112666-30 | Clamshell blister | Klapblister 1 | 90 | 55 | 50 | 40 | 25 | 330 | PET | 1400 | 12 | 16.800 | 13,00 | Gram (g) |
| B112667-30 | Clamshell blister | Klapblister 1 | 90 | 120 | 50 | 100 | 25 | 350 | PET | 600 | 12 | 7.200 | 12,00 | Gram (g) |
| B112667-40 | Clamshell blister | Klapblister 1 | 90 | 120 | 50 | 100 | 35 | 420 | PET | 600 | 12 | 7.200 | 12,00 | Gram (g) |
| B112723-10 | Clamshell blister | Klapblister 1 | 80 | 40 | 55 | 30 | 10 | 310 | PET | 3000 | 12 | 36.000 | 12,00 | Gram (g) |
| B112723-15 | Clamshell blister | Klapblister 1 | 80 | 40 | 55 | 30 | 15 | 310 | PET | 2500 | 12 | 30.000 | 12,00 | Gram (g) |
| B112723-30 | Clamshell blister | Klapblister 1 | 80 | 40 | 55 | 30 | 30 | 420 | PET | 2500 | 12 | 30.000 | 12,00 | Gram (g) |
| B122237-15 | Clamshell blister | Klapblister 1 | 92 | 117 | 65 | 105 | 15 | 310 | PET | 650 | 12 | 7.800 | 13,00 | Gram (g) |
| B122237-25 | Clamshell blister | Klapblister 1 | 92 | 117 | 65 | 105 | 25 | 330 | PET | 650 | 12 | 7.800 | 13,00 | Gram (g) |
| B122237-30 | Clamshell blister | Klapblister 1 | 92 | 117 | 65 | 105 | 30 | 350 | PET | 650 | 12 | 7.800 | 13,00 | Gram (g) |
| B122237-40 | Clamshell blister | Klapblister 1 | 92 | 117 | 65 | 105 | 40 | 420 | PET | 650 | 12 | 7.800 | 13,00 | Gram (g) |
| B122638-40 | Clamshell blister | Klapblister 1 | 195 | 235 | 165 | 190 | 35 | 500 | PET | 100 | 12 | 1.200 | 35,00 | Gram (g) |
| B122638-80 | Clamshell blister | Klapblister 1 | 195 | 235 | 165 | 190 | 75 | 600 | PET | 100 | 12 | 1.200 | 35,00 | Gram (g) |
| B122854-60 | Clamshell blister | Klapblister 1 | 140 | 120 | 100 | 110 | 50 | 600 | PET | 210 | 12 | 2.520 | 14,00 | Gram (g) |

Clappac klapblisters

| | | | | | | | | | | | | | | |
|-----------------|-------------------|---------------|-----|-----|-----|-----|-----|-----|-----|------|----|--------|--------|----------|
| B122854-80 | Clamshell blister | Klapblister 1 | 140 | 120 | 100 | 110 | 70 | 600 | PET | 210 | 12 | 2.520 | 14,00 | Gram (g) |
| B133017 | Clamshell blister | Klapblister 1 | 145 | 80 | 98 | 62 | 15 | 370 | PET | 500 | 15 | 7.500 | 16,00 | Gram (g) |
| B133018 | Clamshell blister | Klapblister 1 | 145 | 80 | 98 | 62 | 30 | 420 | PET | 500 | 15 | 7.500 | 16,00 | Gram (g) |
| B133107-40 | Clamshell blister | Klapblister 1 | 135 | 55 | 105 | 40 | 40 | 420 | PET | 1050 | 12 | 12.600 | 14,00 | Gram (g) |
| CBN10132 | Clamshell blister | Klapblister 1 | 65 | 168 | 40 | 80 | 27 | 400 | PET | 600 | 8 | 4.800 | 8,00 | Gram (g) |
| B133811 | Clamshell blister | Klapblister 1 | 30 | 158 | 28 | 156 | 20 | 280 | PET | 1100 | 5 | 5.500 | 8,00 | Gram (g) |
| CBN10017 5 cavi | Clamshell blister | Klapblister 1 | 110 | 170 | | | | 700 | PET | 470 | 8 | 3.730 | 9,00 | Gram (g) |
| CBN10165 | Clamshell blister | Klapblister 1 | 78 | 120 | 65 | 90 | 40 | 450 | PET | 450 | 28 | 12.600 | 9,40 | Gram (g) |
| CBN10166 | Clamshell blister | Klapblister 1 | 78 | 120 | 65 | 90 | 40 | 400 | PET | 450 | 28 | 12.600 | 9,40 | Gram (g) |
| B113219-100 | Clamshell blister | Klapblister 2 | 175 | 220 | 155 | 177 | 100 | 550 | PET | 120 | 12 | 1.440 | 66,00 | Gram (g) |
| B113219-75 | Clamshell blister | Klapblister 2 | 175 | 220 | 155 | 177 | 75 | 600 | PET | 100 | 15 | 1.500 | 66,00 | Gram (g) |
| B123215-100 | Clamshell blister | Klapblister 2 | 290 | 320 | 250 | 250 | 100 | 500 | PVC | 60 | 12 | 720 | 145,00 | Gram (g) |
| B123217-60 | Clamshell blister | Klapblister 2 | 175 | 436 | 155 | 400 | 60 | 500 | PVC | 65 | 12 | 780 | 100,00 | Gram (g) |
| B123220-90 | Clamshell blister | Klapblister 2 | 175 | 434 | 155 | 400 | 90 | 600 | PVC | 65 | 12 | 780 | 110,00 | Gram (g) |
| B123221-150 | Clamshell blister | Klapblister 2 | 290 | 320 | 247 | 247 | 150 | 600 | PVC | 60 | 12 | 720 | 145,00 | Gram (g) |
| B123223-105 | Clamshell blister | Klapblister 2 | 210 | 655 | 185 | 615 | 105 | 600 | PVC | 35 | 9 | 315 | 188,00 | Gram (g) |
| B123225-60 | Clamshell blister | Klapblister 2 | 210 | 655 | 185 | 615 | 60 | 600 | PVC | 35 | 9 | 315 | 188,00 | Gram (g) |
| B123249-30 | Clamshell blister | Klapblister 2 | 215 | 300 | 190 | 260 | 30 | 350 | PVC | 750 | 2 | 1.500 | 100,00 | Gram (g) |
| B123249-45 | Clamshell blister | Klapblister 2 | 215 | 300 | 190 | 260 | 45 | 450 | PVC | 650 | 1 | 650 | 115,00 | Gram (g) |
| B123249-70 | Clamshell blister | Klapblister 2 | 215 | 300 | 190 | 260 | 70 | 700 | PVC | 600 | 2 | 1.200 | 125,00 | Gram (g) |
| B123373-135 | Clamshell blister | Klapblister 2 | 175 | 434 | 150 | 385 | 135 | 600 | PVC | 60 | 12 | 720 | 110,00 | Gram (g) |
| B132971-20 | Clamshell blister | Klapblister 2 | 130 | 190 | 110 | 150 | 20 | 420 | PET | 300 | 16 | 4.800 | 14,00 | Gram (g) |
| B132971-30 | Clamshell blister | Klapblister 2 | 130 | 190 | 109 | 149 | 30 | 420 | PET | 240 | 12 | 2.880 | 14,00 | Gram (g) |
| B132971-45 | Clamshell blister | Klapblister 2 | 130 | 190 | 107 | 147 | 45 | 580 | PET | 200 | 12 | 2.400 | 14,00 | Gram (g) |
| B133224-30 | Clamshell blister | Klapblister 2 | 140 | 340 | 125 | 297 | 30 | 500 | PVC | 90 | 12 | 1.080 | 40,00 | Gram (g) |
| B133224-60 | Clamshell blister | Klapblister 2 | 145 | 340 | 125 | 297 | 60 | 500 | PET | 90 | 12 | 1.080 | 16,00 | Gram (g) |
| B133226-60 | Clamshell blister | Klapblister 2 | 145 | 170 | 125 | 125 | 60 | 650 | PET | 160 | 9 | 1.440 | 16,00 | Gram (g) |

Clappac klapblisters

| | | | | | | | | | | | | | | |
|---------------|-------------------|----------------|-----|-----|-----|-----|----|-----|-----|------|----|--------|-------|----------|
| B133226-90 | Clamshell blister | Klapblisters 2 | 145 | 170 | 125 | 125 | 90 | 650 | PET | 160 | 9 | 1.440 | 16,00 | Gram (g) |
| B133474-20 | Clamshell blister | Klapblisters 2 | 66 | 120 | 43 | 83 | 20 | 370 | PET | 1250 | 12 | 15.000 | 12,00 | Gram (g) |
| B133474-40 | Clamshell blister | Klapblisters 2 | 66 | 120 | 40 | 80 | 40 | 370 | PET | 1250 | 12 | 15.000 | 12,00 | Gram (g) |
| B133804-20 | Clamshell blister | Klapblisters 2 | 82 | 145 | 61 | 108 | 18 | 250 | PET | 600 | 12 | 7.200 | 18,00 | Gram (g) |
| B133804-30 | Clamshell blister | Klapblisters 2 | 82 | 145 | 61 | 108 | 30 | 800 | PET | 600 | 12 | 7.200 | 18,00 | Gram (g) |
| CBNF/12 | Clamshell blister | Klapblisters 2 | 65 | 145 | 50 | 110 | 12 | 300 | PET | 1150 | 8 | 9.200 | 8,00 | Gram (g) |
| B134105-30 | Clamshell blister | Klapblisters 2 | 140 | 108 | 118 | 70 | 30 | 420 | PET | 400 | 12 | 4.800 | 14,00 | Gram (g) |
| B134105-60 | Clamshell blister | Klapblisters 2 | 140 | 108 | 118 | 70 | 60 | 600 | PET | 500 | 9 | 4.500 | 16,00 | Gram (g) |
| CBNF/20 | Clamshell blister | Klapblisters 2 | 65 | 145 | 50 | 110 | 20 | 400 | PET | 1150 | 8 | 9.200 | 8,00 | Gram (g) |
| CBNF/30 | Clamshell blister | Klapblisters 2 | 65 | 145 | 50 | 110 | 30 | 400 | PET | 1150 | 8 | 9.200 | 8,00 | Gram (g) |
| CBNN01/15 | Clamshell blister | Klapblisters 2 | 67 | 113 | 48 | 75 | 15 | 300 | PET | 1000 | 8 | 8.000 | 7,50 | Gram (g) |
| CBNN01/20 | Clamshell blister | Klapblisters 2 | 67 | 113 | 48 | 75 | 20 | 300 | PET | 1000 | 8 | 8.000 | 7,50 | Gram (g) |
| CBNN01/30 | Clamshell blister | Klapblisters 2 | 67 | 113 | 48 | 75 | 30 | 300 | PET | 1000 | 8 | 8.000 | 7,50 | Gram (g) |
| CBNN01/40 | Clamshell blister | Klapblisters 2 | 67 | 113 | 53 | 75 | 40 | 500 | PET | 1000 | 8 | 8.000 | 7,50 | Gram (g) |
| CBNN01/50 | Clamshell blister | Klapblisters 2 | 67 | 113 | 53 | 75 | 50 | 500 | PET | 1000 | 8 | 8.000 | 7,50 | Gram (g) |
| CBNN02/30 | Clamshell blister | Klapblisters 2 | 85 | 190 | 62 | 150 | 30 | 400 | PET | 800 | 8 | 6.400 | 12,00 | Gram (g) |
| CBNN03/30 | Clamshell blister | Klapblisters 2 | 98 | 200 | 72 | 160 | 30 | 400 | PET | 500 | 8 | 4.000 | 15,00 | Gram (g) |
| CBNN03/45 | Clamshell blister | Klapblisters 2 | 98 | 200 | 70 | 160 | 45 | 400 | PET | 500 | 8 | 4.000 | 15,00 | Gram (g) |
| CBNN04/20 | Clamshell blister | Klapblisters 2 | 88 | 165 | 66 | 130 | 20 | 400 | PET | 600 | 8 | 4.800 | 10,00 | Gram (g) |
| CBNN04/6 cavi | Clamshell blister | Klapblisters 2 | 88 | 165 | | | | | PET | 600 | 8 | 4.800 | | Gram (g) |
| CBNN100/12 | Clamshell blister | Klapblisters 2 | 67 | 80 | 50 | 50 | 12 | 300 | PET | 1800 | 8 | 14.400 | 6,00 | Gram (g) |
| CBNN100/20 | Clamshell blister | Klapblisters 2 | 67 | 80 | 50 | 50 | 20 | 300 | PET | 1800 | 8 | 14.400 | 6,00 | Gram (g) |
| CBNN100/30 | Clamshell blister | Klapblisters 2 | 67 | 80 | 50 | 50 | 30 | 400 | PET | 1800 | 8 | 14.400 | 6,00 | Gram (g) |
| CBNN100/40 | Clamshell blister | Klapblisters 2 | 67 | 80 | 50 | 50 | 40 | 400 | PET | 1800 | 8 | 14.400 | 6,00 | Gram (g) |
| CBNN110/10 | Clamshell blister | Klapblisters 2 | 340 | 100 | 310 | 75 | 12 | 300 | PET | 300 | 8 | 2.400 | 25,00 | Gram (g) |
| CBNN110/30 | Clamshell blister | Klapblisters 2 | 340 | 100 | 310 | 75 | 30 | 400 | PET | 300 | 8 | 2.400 | 25,00 | Gram (g) |
| CBNN110/50 | Clamshell blister | Klapblisters 2 | 340 | 100 | 310 | 75 | 50 | 500 | PET | 300 | 8 | 2.400 | 25,00 | Gram (g) |
| CBNN120/12 | Clamshell blister | Klapblisters 2 | 68 | 175 | 50 | 150 | 12 | 300 | PET | 700 | 8 | 5.600 | 8,20 | Gram (g) |
| CBNN120/20 | Clamshell blister | Klapblisters 2 | 68 | 175 | 50 | 150 | 20 | 300 | PET | 700 | 8 | 5.600 | 8,20 | Gram (g) |
| CBNN120/40 | Clamshell blister | Klapblisters 2 | 68 | 175 | 48 | 141 | 40 | 400 | PET | 700 | 8 | 5.600 | 8,20 | Gram (g) |
| CBNN130/10 | Clamshell blister | Klapblisters 2 | 73 | 130 | 65 | 95 | 10 | 400 | PET | 930 | 8 | 7.440 | 9,40 | Gram (g) |
| CBNN140/30 | Clamshell blister | Klapblisters 2 | 92 | 140 | 75 | 105 | 30 | 400 | PET | 1000 | 8 | 8.000 | 12,00 | Gram (g) |
| CBNN140/40 | Clamshell blister | Klapblisters 2 | 92 | 140 | 75 | 105 | 40 | 500 | PET | 1000 | 8 | 8.000 | 12,00 | Gram (g) |
| CBNN140/60 | Clamshell blister | Klapblisters 2 | 92 | 140 | 72 | 100 | 60 | 500 | PET | 800 | 8 | 6.400 | 12,00 | Gram (g) |
| CBNN150/30 | Clamshell blister | Klapblisters 2 | 94 | 260 | 75 | 220 | 30 | 400 | PET | 700 | 8 | 5.600 | 16,00 | Gram (g) |

Clappac klapblisters

| | | | | | | | | | | | | | | |
|------------|-------------------|----------------|-----|-----|-----|-----|----|-----|-----|-----|----|-------|-------|----------|
| CBNN150/40 | Clamshell blister | Klapblisters 2 | 94 | 260 | 72 | 212 | 40 | 400 | PET | 700 | 8 | 5.600 | 16,00 | Gram (g) |
| CBNN160/20 | Clamshell blister | Klapblisters 2 | 86 | 129 | 62 | 86 | 20 | 300 | PET | 900 | 8 | 7.200 | 10,00 | Gram (g) |
| CBNN160/35 | Clamshell blister | Klapblisters 2 | 125 | 86 | 86 | 63 | 35 | 400 | PET | 850 | 8 | 6.800 | 12,00 | Gram (g) |
| CBNN170/30 | Clamshell blister | Klapblisters 2 | 102 | 120 | 80 | 80 | 30 | 400 | PET | 740 | 8 | 5.920 | 9,00 | Gram (g) |
| CBNN170/50 | Clamshell blister | Klapblisters 2 | 102 | 120 | 78 | 78 | 50 | 500 | PET | 740 | 8 | 5.920 | 9,00 | Gram (g) |
| CBNN180/30 | Clamshell blister | Klapblisters 2 | 110 | 190 | 90 | 150 | 30 | 400 | PET | 500 | 8 | 4.000 | 9,00 | Gram (g) |
| CBNN180/50 | Clamshell blister | Klapblisters 2 | 110 | 190 | 90 | 150 | 50 | 500 | PET | 500 | 8 | 4.000 | 9,00 | Gram (g) |
| CBNN180/60 | Clamshell blister | Klapblisters 2 | 110 | 190 | 90 | 150 | 60 | 500 | PET | 500 | 8 | 4.000 | 9,00 | Gram (g) |
| CBNN190/35 | Clamshell blister | Klapblisters 2 | 110 | 190 | 115 | 155 | 35 | 400 | PET | 180 | 16 | 2.880 | 9,00 | Gram (g) |
| CBNN190/60 | Clamshell blister | Klapblisters 2 | 110 | 190 | 112 | 152 | 60 | 500 | PET | 360 | 8 | 2.880 | 9,00 | Gram (g) |
| CBNN230/40 | Clamshell blister | Klapblisters 2 | 140 | 240 | 120 | 200 | 40 | 400 | PET | 250 | 8 | 2.000 | 17,00 | Gram (g) |
| CBNN230/60 | Clamshell blister | Klapblisters 2 | 140 | 240 | 115 | 195 | 60 | 400 | PET | 250 | 8 | 2.000 | 21,00 | Gram (g) |

| | | | | | | | | | | | | | | |
|----------------|-------------------|----------------|-----|-----|-----|-----|----|-----|-----|-----|----|-------|-------|----------|
| CBNN240/30 | Clamshell blister | Klapblisters 2 | 142 | 120 | 125 | 90 | 30 | 400 | PET | 600 | 8 | 4.800 | 19,00 | Gram (g) |
| CBNN240/40 | Clamshell blister | Klapblisters 2 | 142 | 120 | 121 | 88 | 40 | 400 | PET | 600 | 8 | 4.800 | 19,00 | Gram (g) |
| CBNN250/30 | Clamshell blister | Klapblisters 2 | 148 | 170 | 125 | 125 | 30 | 400 | PET | 400 | 8 | 3.200 | 19,00 | Gram (g) |
| CBNN250/40 | Clamshell blister | Klapblisters 2 | 148 | 170 | 125 | 125 | 40 | 400 | PET | 400 | 8 | 3.200 | 19,00 | Gram (g) |
| CBNN250/50 | Clamshell blister | Klapblisters 2 | 148 | 170 | 125 | 125 | 50 | 500 | PET | 400 | 8 | 3.200 | 19,00 | Gram (g) |
| CBNN250/90 | Clamshell blister | Klapblisters 2 | 148 | 170 | 121 | 121 | 90 | 700 | PET | 300 | 8 | 2.400 | 18,00 | Gram (g) |
| CBNN270/20 | Clamshell blister | Klapblisters 2 | 190 | 260 | 165 | 215 | 20 | 300 | PET | 160 | 8 | 1.280 | 32,00 | Gram (g) |
| CBNN270/50 | Clamshell blister | Klapblisters 2 | 190 | 260 | 165 | 215 | 50 | 500 | PET | 160 | 8 | 1.280 | 32,00 | Gram (g) |
| CBNN270/60 | Clamshell blister | Klapblisters 2 | 190 | 260 | 165 | 215 | 60 | 600 | PET | 160 | 8 | 1.280 | 32,00 | Gram (g) |
| CBNN280/30 | Clamshell blister | Klapblisters 2 | 188 | 115 | 167 | 80 | 30 | 400 | PET | 400 | 8 | 3.200 | 32,00 | Gram (g) |
| CBNN280/40 | Clamshell blister | Klapblisters 2 | 188 | 115 | 167 | 78 | 40 | 400 | PET | 400 | 8 | 3.200 | 32,00 | Gram (g) |
| CBNN280/70 | Clamshell blister | Klapblisters 2 | 188 | 115 | 165 | 75 | 70 | 600 | PET | 360 | 8 | 2.880 | 32,00 | Gram (g) |
| CBNN320/40 | Clamshell blister | Klapblisters 2 | 235 | 290 | 200 | 240 | 40 | 400 | PET | 130 | 10 | 1.300 | 25,00 | Gram (g) |
| CBNN320/80 | Clamshell blister | Klapblisters 2 | 235 | 290 | 190 | 230 | 80 | 600 | PET | 100 | 10 | 1.000 | 25,00 | Gram (g) |
| B133218-95 | Clamshell blister | Klapblisters 2 | 165 | 190 | 145 | 145 | 95 | 500 | PET | 130 | 12 | 1.560 | 16,00 | Gram (g) |
| B133425 | Clamshell blister | Klapblisters 2 | 165 | 165 | 125 | 145 | 10 | 350 | PET | 200 | 12 | 2.400 | 16,00 | Gram (g) |
| CBNN260/20 | Clamshell blister | Klapblisters 2 | 165 | 190 | 145 | 148 | 20 | 300 | PET | 300 | 8 | 2.400 | 30,00 | Gram (g) |
| CBNN260/30 | Clamshell blister | Klapblisters 2 | 165 | 190 | 145 | 148 | 30 | 400 | PET | 300 | 8 | 2.400 | 30,00 | Gram (g) |
| CBN10023 | Clamshell blister | Klapblisters 3 | 65 | 168 | 40 | 60 | 27 | 400 | PET | 600 | 16 | 9.600 | 8,00 | Gram (g) |
| B133716-30 | Clamshell blister | Klapblisters 3 | 105 | 145 | 80 | 80 | 30 | 370 | PET | 480 | 12 | 5.760 | 18,00 | Gram (g) |
| B133716-504025 | Clamshell blister | Klapblisters 3 | 105 | 145 | 40 | 50 | 25 | 370 | PET | 480 | 12 | 5.760 | 18,00 | Gram (g) |
| B133716-656520 | Clamshell blister | Klapblisters 3 | 105 | 145 | 65 | 65 | 20 | 370 | PET | 480 | 12 | 5.760 | 18,00 | Gram (g) |

Clappac klapblisters

| | | | | | | | | | | | | | | |
|--------------------|-------------------|----------------|-----|-----|-----|-----|----|-----|-----|------|----|--------|-------|----------|
| B133716-656530 | Clamshell blister | Klapblisters 3 | 105 | 145 | 65 | 65 | 30 | 370 | PET | 480 | 12 | 5.760 | 18,00 | Gram (g) |
| B133716-808035 | Clamshell blister | Klapblisters 3 | 105 | 145 | 80 | 80 | 35 | 370 | PET | 480 | 12 | 5.760 | 18,00 | Gram (g) |
| B133717-1558035 | Clamshell blister | Klapblisters 3 | 108 | 243 | 80 | 155 | 35 | 400 | PET | 300 | 12 | 3.600 | 19,00 | Gram (g) |
| B133717-1558070 | Clamshell blister | Klapblisters 3 | 108 | 243 | 80 | 155 | 70 | 550 | PET | 280 | 11 | 3.120 | 19,00 | Gram (g) |
| B133717-207 | Clamshell blister | Klapblisters 3 | 108 | 243 | 78 | 203 | 29 | 400 | PET | 300 | 12 | 3.600 | 19,00 | Gram (g) |
| B133717-30 | Clamshell blister | Klapblisters 3 | 108 | 243 | 80 | 180 | 30 | 400 | PET | 300 | 12 | 3.600 | 20,00 | Gram (g) |
| B133718-35 | Clamshell blister | Klapblisters 3 | 190 | 204 | 155 | 120 | 35 | 420 | PET | 200 | 10 | 2.000 | 18,00 | Gram (g) |
| B133718-55 | Clamshell blister | Klapblisters 3 | 190 | 204 | 155 | 120 | 55 | 550 | PET | 150 | 12 | 1.800 | 18,00 | Gram (g) |
| B133725-1338035 | Clamshell blister | Klapblisters 3 | 110 | 220 | 80 | 133 | 35 | 400 | PET | 280 | 12 | 3.360 | 20,00 | Gram (g) |
| B133725-1338045 | Clamshell blister | Klapblisters 3 | 110 | 220 | 80 | 133 | 45 | 450 | PET | 280 | 12 | 3.360 | 20,00 | Gram (g) |
| B133725-1659025 | Clamshell blister | Klapblisters 3 | 110 | 220 | 90 | 165 | 25 | 400 | PET | 330 | 16 | 5.280 | 20,00 | Gram (g) |
| B133725-23 | Clamshell blister | Klapblisters 3 | 110 | 220 | 80 | 155 | 25 | 400 | PET | 280 | 12 | 3.360 | 20,00 | Gram (g) |
| B133725-35 | Clamshell blister | Klapblisters 3 | 110 | 220 | 80 | 155 | 35 | 400 | PET | 280 | 12 | 3.360 | 20,00 | Gram (g) |
| B133725-70 | Clamshell blister | Klapblisters 3 | 110 | 220 | 80 | 155 | 70 | 550 | PET | 240 | 12 | 2.880 | 20,00 | Gram (g) |
| B133760-14414437 | Clamshell blister | Klapblisters 3 | 107 | 180 | 144 | 144 | 35 | 370 | PET | 200 | 12 | 2.400 | 26,00 | Gram (g) |
| B133772 - 13013074 | Clamshell blister | Klapblisters 3 | 155 | 190 | 130 | 130 | 74 | 800 | PET | 500 | 9 | 4.500 | 16,00 | Gram (g) |
| B133804-42 | Clamshell blister | Klapblisters 3 | 82 | 145 | 61 | 108 | 42 | 800 | PET | 400 | 16 | 6.400 | 19,00 | Gram (g) |
| B134013-23 | Clamshell blister | Klapblisters 3 | 110 | 220 | 80 | 80 | 25 | 400 | PET | 280 | 12 | 3.360 | 20,00 | Gram (g) |
| B134176-1104520 | Clamshell blister | Klapblisters 3 | 90 | 180 | 45 | 110 | 20 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1206515 | Clamshell blister | Klapblisters 3 | 90 | 180 | 65 | 120 | 15 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1206520 | Clamshell blister | Klapblisters 3 | 90 | 180 | 65 | 120 | 20 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1206535 | Clamshell blister | Klapblisters 3 | 90 | 180 | 65 | 120 | 35 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1307515 | Clamshell blister | Klapblisters 3 | 90 | 180 | 75 | 130 | 15 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1307525 | Clamshell blister | Klapblisters 3 | 90 | 180 | 75 | 130 | 25 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1307535 | Clamshell blister | Klapblisters 3 | 90 | 180 | 75 | 130 | 35 | 400 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1404520 | Clamshell blister | Klapblisters 3 | 90 | 180 | 45 | 140 | 20 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1457715 | Clamshell blister | Klapblisters 3 | 90 | 180 | 77 | 145 | 15 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-1457725 | Clamshell blister | Klapblisters 3 | 90 | 180 | 77 | 145 | 25 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-806515 | Clamshell blister | Klapblisters 3 | 90 | 180 | 80 | 65 | 15 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134176-967530 | Clamshell blister | Klapblisters 3 | 90 | 180 | 96 | 75 | 30 | 350 | PET | 250 | 22 | 5.500 | 14,00 | Gram (g) |
| B134183-405040 | Clamshell blister | Klapblisters 3 | 75 | 135 | 40 | 50 | 40 | 400 | PET | 900 | 14 | 12.600 | 12,00 | Gram (g) |
| B134183-504015 | Clamshell blister | Klapblisters 3 | 75 | 135 | 40 | 50 | 15 | 330 | PET | 1000 | 14 | 14.000 | 12,00 | Gram (g) |
| B134183-504025 | Clamshell blister | Klapblisters 3 | 75 | 135 | 40 | 50 | 25 | 330 | PET | 900 | 14 | 12.600 | 12,00 | Gram (g) |
| B134183-705520 | Clamshell blister | Klapblisters 3 | 75 | 135 | 55 | 70 | 20 | 330 | PET | 1000 | 12 | 12.000 | 12,00 | Gram (g) |

Clappac klapblisters

| | | | | | | | | | | | | | | |
|-----------------|-------------------|----------------|-----|-----|----|-----|----|-----|-----|------|----|--------|-------|----------|
| B134183-856025 | Clamshell blister | Klapblisters 3 | 75 | 135 | 60 | 85 | 25 | 400 | PET | 900 | 14 | 12.600 | 12,00 | Gram (g) |
| B134183-856040 | Clamshell blister | Klapblisters 3 | 75 | 135 | 60 | 85 | 40 | 400 | PET | 800 | 12 | 9.600 | 12,00 | Gram (g) |
| B134222-504015 | Clamshell blister | Klapblisters 3 | 75 | 100 | 50 | 40 | 15 | 350 | PET | 1000 | 12 | 12.000 | 12,00 | Gram (g) |
| B134230-1309030 | Clamshell blister | Klapblisters 3 | 107 | 180 | 90 | 130 | 30 | 420 | PET | 200 | 18 | 3.600 | 19,00 | Gram (g) |
| B134230-1309040 | Clamshell blister | Klapblisters 3 | 107 | 180 | 90 | 130 | 40 | 420 | PET | 520 | 9 | 4.680 | 19,00 | Gram (g) |
| B134230-1459030 | Clamshell blister | Klapblisters 3 | 180 | 107 | 90 | 145 | 30 | 420 | PET | 200 | 18 | 3.600 | 18,00 | Gram (g) |
| B134281-784510 | Clamshell blister | Klapblisters 3 | 60 | 105 | 45 | 78 | 10 | 400 | PET | 1400 | 12 | 16.800 | 12,00 | Gram (g) |
| B134281-784515 | Clamshell blister | Klapblisters 3 | 60 | 105 | 45 | 78 | 15 | 400 | PET | 1400 | 12 | 16.800 | 12,00 | Gram (g) |
| B134281-784525 | Clamshell blister | Klapblisters 3 | 60 | 105 | 45 | 75 | 25 | 400 | PET | 1400 | 12 | 16.800 | 12,00 | Gram (g) |
| B134309-1104025 | Clamshell blister | Klapblisters 3 | 60 | 150 | 40 | 110 | 25 | 420 | PET | 900 | 12 | 10.800 | 12,00 | Gram (g) |
| B134309-1104030 | Clamshell blister | Klapblisters 3 | 60 | 150 | 40 | 110 | 30 | 420 | PET | 900 | 12 | 10.800 | 12,00 | Gram (g) |
| CBN10005 | Clamshell blister | Klapblisters 3 | 95 | 187 | 70 | 103 | 25 | 400 | PET | 600 | 8 | 4.800 | 16,00 | Gram (g) |
| CBN10006 | Clamshell blister | Klapblisters 3 | 95 | 187 | 55 | 82 | 36 | 400 | PET | 600 | 8 | 4.800 | 11,00 | Gram (g) |
| CBN10007 | Clamshell blister | Klapblisters 3 | 95 | 187 | 60 | 80 | 50 | 400 | PET | 600 | 8 | 4.800 | 11,00 | Gram (g) |
| CBN10066 USB | Clamshell blister | Klapblisters 3 | 95 | 148 | 18 | 70 | 10 | 400 | PET | 480 | 16 | 7.680 | 16,00 | Gram (g) |
| CBN10084 | Clamshell blister | Klapblisters 3 | 100 | 168 | 73 | 73 | 20 | 400 | PET | 400 | 16 | 6.400 | 5,00 | Gram (g) |
| CBN10092 | Clamshell blister | Klapblisters 3 | 100 | 168 | 73 | 33 | 20 | 400 | PET | 400 | 16 | 6.400 | 5,00 | Gram (g) |
| CBN10134 | Clamshell blister | Klapblisters 3 | 98 | 190 | 78 | 152 | 20 | 400 | PET | 600 | 8 | 4.800 | 15,00 | Gram (g) |
| CBN10143 | Clamshell blister | Klapblisters 3 | 110 | 202 | 90 | 165 | 20 | 400 | PET | 600 | 8 | 4.800 | 9,00 | Gram (g) |
| CBN4097 | Clamshell blister | Klapblisters 3 | 95 | 187 | 70 | 97 | 50 | 400 | PET | 600 | 8 | 4.800 | 11,00 | Gram (g) |
| CBN4098 | Clamshell blister | Klapblisters 3 | 145 | 173 | 85 | 108 | 43 | 500 | PET | 550 | 8 | 4.400 | 19,00 | Gram (g) |
| CBNN01/1 | Clamshell blister | Klapblisters 3 | 67 | 113 | 45 | 35 | 20 | 300 | PET | 1000 | 8 | 8.000 | 7,50 | Gram (g) |
| CBNN01/10 | Clamshell blister | Klapblisters 3 | 67 | 113 | 45 | 35 | 10 | 300 | PET | 1000 | 8 | 8.000 | 7,50 | Gram (g) |
| B134177-1109030 | Clamshell blister | Klapblisters 3 | 115 | 155 | 90 | 110 | 30 | 350 | PET | 280 | 22 | 6.160 | 20,00 | Gram (g) |